


| | | | |
|--|--|----------------------------|--|
|  | Public Health Reference Laboratory (PHRL) | |  |
| | Patient Investigating Form for Mumps | | |
| Page 1 of 8 | Doc# KMU/PHRL/SOP/24 | Effective Date: 07-06-2022 | Version: 01 |

| | | |
|-----|---------------------|--|
| 01. | Patient Name | |
| 02. | Father's Name | |
| 03. | Age | |
| 04. | CNIC | |
| 05. | Contact Number | |
| 06. | Referring Physician | |
| 07. | Address | |

Referred test: Anti Mumps (ELISA)

History:

Comments:

SOPs to be followed:

1. Ship specimens for testing as soon as possible.
2. Samples must be collected in gel/serum separator tube (tiger top tube)
3. Clearly label the tube with a permanent marker.
4. Complete the patient investigation form.
5. Pack the tube in proper packaging.
6. Cover spill completely with paper towels and pour disinfectant solution (0.5% sodium hypochlorite) over paper towels. Allow paper towels soaked with disinfectant to stand for 30 minutes. Clean up the spill. Rinse area with 70% ethanol or water. Place spill clean-up materials in a biomedical waste container (Yellow bag). Wash hands with soap and water.
7. Include a frozen cold pack with the specimen.
8. Write the address on the packaging.
9. Surgical mask and gloves are recommended for all procedure.


⌘-----

ELISA Section

Immunology and Serology Lab

PHRL, KMU phase 5, Peshawar.

Contact number:

| | | | |
|--|--|----------------------------|--|
|  | Public Health Reference Laboratory (PHRL) | |  |
| | Mumps Description | | |
| Page 2 of 8 | Doc# KMU/PHRL/SOP/24 | Effective Date: 07-06-2022 | Version: 01 |

Mumps

Mumps is a virus-borne infectious illness. It usually begins with a fever, headache, muscular pains, exhaustion, and a loss of appetite that lasts a few days. The salivary glands of most people enlarge (often referred to as parotitis when the parotid gland, located in front and below the ear, swells). The plump cheeks and sensitive, swelling jaw are the result of this.

Signs & Symptoms of Mumps

Mumps is most well-known for causing puffy cheeks and a painful, swollen jaw. This is caused by swelling salivary glands under one or both ears, a condition known as parotitis.

Other signs and symptoms that may appear a few days before the onset of parotitis include:

- ❖ Fever
- ❖ Headache
- ❖ Muscle aches
- ❖ Tiredness
- ❖ Loss of appetite

Symptoms typically appear 16-18 days after infection, but this period can range from 12–25 days after infection.



Mumps can cause moderate symptoms (like a cold) or no symptoms at all, and some individuals are unaware they have the condition.

Mumps can sometimes lead to more serious issues.

Mumps usually clear up with in two weeks for the majority of people.

Transmission of Mumps

Mumps is a virus-borne infectious illness. It spreads when saliva or respiratory droplets from the lips, nose, or throat come into touch with it. An infected person can spread the virus by

- ❖ coughing, sneezing, or talking
- ❖ sharing items that may have saliva on them, such as water bottles or cups
- ❖ participating in close-contact activities with others, such as playing sports, dancing, or kissing

Mumps can be disseminated from a few days before the salivary glands swell to up to five days after the swelling starts. During this period, somebody with mumps should avoid interacting with others. Stay home from school and avoid social gatherings.

Use of Serology to Aid in the Diagnosis of Mumps Infection


To confirm mumps infection, a buccal swab specimen taken 3 days after the beginning of parotitis is preferable. Detection of mumps immunoglobulin M (IgM) can aid in the diagnosis of mumps although a positive IgM result determines a probable rather than confirmed case, based on the Council of State and Territorial Epidemiologists mumps case definition. A significant rise in immunoglobulin G (IgG) antibody titer in paired acute- and convalescent-phase serum specimens, and IgG seroconversion, can also aid in the diagnosis of mumps infection but are not recommended because they are rarely observed.

The ability to detect IgM in people infected with mumps is influenced by vaccination status and the timing of test collection. IgM detection is highest in unprotected people, middle in one-dose vaccination receivers, and lowest in two-dose vaccine recipients in general. A positive IgM result is more common in samples obtained more than 3 days after the beginning of parotitis.

The sensitivity and specificity of IgM test techniques and kits varies greatly, with certain indirect EIA and immunofluorescent tests identifying as little as 12–15 percent of confirmed mumps infections.

The IgM capture ELISA is the most sensitive serologic approach, identifying 46 percent to 71% of rRT PCR confirmed infections, although it is only available in limited quantities.

It has been widely known that mumps IgM cannot be detected in previously vaccinated people who get sick with the disease. Regardless of when the material is collected, those who have had mumps immunization may not have detectable mumps IgM antibodies. Because IgG test results are usually positive and raised at the time of the initial blood collection, detecting a 4-fold spike or seroconversion is improbable and hence not advised.

| | | | |
|--|--|----------------------------|--|
|  | Public Health Reference Laboratory (PHRL) | |  |
| | SOP–Mumps Sample Collection | | |
| Page 4 of 8 | Doc# KMU/PHRL/SOP/24 | Effective Date: 07-06-2022 | Version: 01 |

1.0 Purpose

To maintain patient safety by providing consistent method for correct identification, sample collection, transportation, receiving and/or rejection of specimens of suspected patients for Mumps.

2.0 Scope

This policy applies to all staff and areas of Clinical Laboratory where samples are collected, received, registered, processed, archive and disposal by maintaining standards.

3.0 Responsibility

- 3.1 All staff receiving samples must comply with the procedure in acceptance and rejection of samples.
- 3.2 Lab Head and Managers/Assistant Managers/Coordinator are responsible for implementation and monitoring of policies/procedures.
- 3.3 Managers/Assistant Managers/Coordinators are also responsible for providing awareness and understanding of policy.
- 3.4 Senior Coordinator, QA is responsible for revision and updating of policy.
- 3.5 Laboratory Director and Senior Administrator are responsible for the approval.

4.0 Procedure

- 4.1 Proper Personal Protective Equipment should be used while handling biological specimen.
- 4.2 Specimens are transported to receiving area in individual plastic bags (zip lock bags) in a labeled transportation box accompanied by requisition slip.
- 4.3 Upon arrival following steps must be taken:
 - 4.3.1 Remove the sample from the rack or bag as per required safety precautions.
 - 4.3.2 Check that the specimen tube/container is firmly capped and with no leakage.
 - 4.3.3 Store samples at room temperature for up to 24 hours.
 - 4.3.4 If you need to store the sample for one week before shipping, store between 0-5° Celsius.
 - 4.3.5 Avoid freezing whole blood samples.
 - 4.3.6 Read the patient's full name and L# from the requisition slip.
 - 4.3.7 Check the sample is labeled with same name and identification number.

- 4.3.8** If the sample is unlabeled or poorly labeled treat it as rejected sample (see rejection criteria Policy Document No. KMU/PHRL/SOP-04)
- 4.3.9** All patient information provided on requisition slip should be clear, unambiguous and allow for unequivocal identification of the patient, the requesting physician, patient location, clinical details and investigations where required.
- 4.4** Technologist must read the requisition slip for any additional comments/remarks for example STAT test (STAT laboratory tests and services are those that are needed immediately to manage medical emergencies) and request for special processing.
- 4.5** Samples received from outreach collection units at outreach Lab are check for the patient's full name and MR# or Lab# from the requisition slip and check the sample is labeled with same name and identification number then distributed to respective sub-section/ benches.
- 4.6** Specimens are then sorted according to laboratory sub-section on priority basis. If test request received marked with "STAT" sample must be delivered to sub-section immediately with special attention and test should be carried out within specified time of STAT tests.
- 4.7** If necessary samples are spun and aliquot prepared.
- 4.8** Depending upon the tests requested, all samples are may be stored at room temperature (15°-30°C) for up to two hours prior to testing or transferred immediately to clod chain storage of other diseases testing lab
- 4.9** If the sample is satisfactory and need to be shared with any other sub section/benches of clinical laboratory, make aliquots with proper identification and it should be sent immediately to relevant sub section/benches along with requisition slip.
- 4.10** If sample is accepted in the outreach laboratory or Collection center, it is the responsibility of outreach laboratory staff to protect the sample according to the outreach laboratory retention policy.
- 4.11** Receiving bench technologists are responsible for taking out list of samples from LIMS/IPMS that is collected and not received in laboratory; after defined period of time it should be informed to relevant area supervisor (Reception or Laboratory

Outreach Support Centre) immediately for retrieval of specimen to avoid delay in processing and reporting.

4.12 All unsatisfactory samples entered in the computer record are immediately cancelled stating the reason for rejection. All such information is documented in the comments section of the LIMS (Laboratory Information Management System).

4.13 Unsatisfactory sample will not be returned back to in-patient area and is stored according to the sample retention policy of Outreach Lab.


5.0 General Notes for Handling Reagents and/or Samples

- 5.1** Handle all specimens as if they contain infectious agents
- 5.2** Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper disposal of specimens.
- 5.3** The work area/benches are disinfected before and after the test(s) performed to avoid any bio-hazards risk.
- 5.4** During processing samples, staff must wear proper PPE (Personal protective equipment).
- 5.5** Mouth pipetting of any chemical or sample material is completely prohibited.
- 5.6** Precaution should be taken to avoid punctures, cuts and any open skin wounds.
- 5.7** Adhere and follow safety guidelines when handling chemicals and infectious sample materials.
- 5.8** Contain waste appropriately for disposal.
- 5.9** Gloves should be removed in such a way that the skin does not come into contact with the external surface of the glove.
- 5.10** Do not handle taps, phone, switches and door handle directly, while working with infectious sample.
- 5.11** Wash hand thoroughly after processing samples.

6.0 REFERENCE

7.0 APPROVAL

| | |
|--------------------|------------------|
| Prepared by | |
| Name & Designation | Signature & Date |
| Reviewed by | |
| Name & Designation | Signature & Date |
| Name & Designation | Signature & Date |
| Approved by | |
| Name & Designation | Signature & Date |

| | | | |
|--|---|----------------------------|--|
|  | Public Health Reference Laboratory (PHRL) | |  |
| | SOP –Mumps Sample Packing and Transportation to PHRL-KMU. | | |
| Page 7 of 2 | Doc# KMU/PHRL/SOPs/24 | Effective Date: 07-06-2022 | Version: 01 |

Public Health Reference Laboratory (PHRL) is committed to provide quality services for disease surveillance through Integrated Disease Surveillance and Response (IDSR) in Khyber Pakhtunkhwa province.

8.0 Purpose

To describe the laboratory procedures for Mumps Sample Packing and Transportation to PHRL, KMU.

9.0 Scope

These Standard Operating Procedures provide a template agreed by PHRLs that is essential for guiding responses to infectious disease outbreaks by public health practitioners, epidemiologists and clinicians.

10.0 Responsibility

- 10.1** Doing: Lab in charge/Team lead/Technicians.
- 10.2** Checking/Reporting: Consultants/Reporting group.
- 10.3** All trained personnel are responsible to handle and test samples suspected of containing Mumps.
- 10.4** Laboratory Director/Senior Administrator are responsible for the approval, implementation, validation, maintenance and review of this procedure.

11.0 Procedure

11.1 Specimen Packing and Transportation

- 11.1.1** Samples container (blood/sera tubes) must be properly labelled with double identification marks i.e. NAME & Sample ID perfectly reflecting the excel sheet attached.
- 11.1.2** Absorbent materials should be in sufficient quantity to absorb the entire liquid content if the primary container (s) leak.
- 11.1.3** Leak proof secondary container, cushioning material e.g. bubble wrap, tape to seal the outer package (if required).

- 11.1.4 Excel sheets of the samples must be double checked for any duplication, in case of similarity in names.
- 11.1.5 Laboratory form or letter describing the main epidemiological and clinical findings and the lab tests that are required
- 11.1.6 Make sure that samples must be sorted in small groups of 05-10/ bag, with individual requisition/Referral form for all samples attached/inserted in the pocket.
- 11.1.7 After collection and packaging the sample, it must be immediately sent to PHRL-KMU.
- 11.1.8 Receiver name, address and telephone number.

11.2 Specimen Rejection Criteria

- 11.2.1 If specimen does not meet acceptable criteria it may be rejected by the laboratory. The following are criteria used for the possible rejection of specimens:
- 11.2.2 Specimen improperly labeled as to the patient identity (as per above criteria)
- 11.2.3 Patient identification mismatched between specimen and requisition form/ excel sheet data.
- 11.2.4 Improper collection container i.e., leakage or broken sample tube or container.
- 11.2.5 Sample collection date is older than 7 days from the date of receiving.
- 11.2.6 Improperly transported (i.e., not on required temperature)
- 11.2.7 Specimen contaminated with biological hazardous material.

General Notes

- All Labs can coordinate and follow-up the results on 091-9219651 and each Lab are expected to maintain complete record of samples sending, results and communication with PHRL focal person.
- For any further assistance Mr. Abdul Qadoos khan (Lab Technologist) may be contacted through email:aqkhan.kmu@gmail.com